



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Dental applicators KERR B 298 (black) (New)

New,
U.S.A. production,
Kerr applicators are used to apply binding systems used in dentistry, e.g:
Die-cutters,
Binding materials,
Filling primers,
Fluoride varnish,
Liquid to stop the bleeding,
Liquids for the determination of caries outbreaks,
Disposable,
They are made of high quality flexible material,
The applicator tip is made of fibre,
Length of applicator: 8.5 cm,
Package dimensions: 10,4 x 2,2 x 2,2 cm,
The package contains 50 applicators,

