

PRAISTON Sp. z o.o. ul. Górowska 32 64-100 Leszno Polska Biuro: +48 65 527 01 67 Serwis:+48 782 844 000 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży Krzysztof Wybieralski +48 882 762 006

kw@praiston.pl



RUSCH LIFE-AIR 1000 Patient warmer (Reconditioned)

- Used,
- · Technical condition: very good,
- · Visual condition: good,
- · Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- LIFE AIR 1000 hypothermic therapy system was developed to reduce the effects of hypothermia, which is a common problem during and after surgery,
- Manual control of the temperature of air flowing around the patient with the heater control switch in a selectable range:
 - 。27° C
 - ° 32° C +/- 3° C,
 - ° 36° C +/- 3° C,
 - 43° C +/- 3° C,
- · Equipped with:
 - Control panel in English,
 - · Test lamp built-in automatic test sequence to check the correct operation of heating elements,
 - Warning light visual and audible alarm to protect the system from overheating,
- Temperature Bar Graph visually displays the "range" of the temperature of the air administered to the patient. Provides the operator with a quick visual check of the temperature range,
 - Power Indicator,
- · Low weight for easy transport and preparation,
- · Mobile on double casters,
- Includes:
 - User manual in English(PDF),
- · Dimensions:
 - Heater: 29 x 30.5 x 38/58 cm,
 - Overall: 53 x 53 x 75 cm,
- Weight: 15.5 kg,
- · Valid technical passport issued,
- Warranty:
 - · 6 months for domestic market Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU. WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски Мы говорим по украински



www.praiston.pl