



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **EYE-LIGHT Dry eye treatment device (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Italy, (2018),
- The unit worked only 107 hours,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Exclusive and patented technologies to optimize the use of light and heat energy,
- OPE is a patented broad-spectrum light technology that excludes all emissions irrelevant to treatment. It is based on a highly customizable technology capable of optimizing tissue stimulation and regeneration,
- Thanks to a patented and unique cooling module, no gel is used, making the treatment simpler, faster and less expensive,
- Highly intuitive software allows the operator to work optimally and safely from the first treatment, with automatic adjustment of the power and duration of each light pulse.
- LIGHT MODULATION - LLLT / LOW LEVEL LIGHT THERAPY:
  - LIGHT MODULATION® is a patented photobiomodulation (PBM) technology, a powerful metabolic enhancer that stimulates ATP (adenosine triphosphate) production to increase cellular action and accentuate cellular activity,
  - The LED emitting matrix releases endogenous eyelid heat, directly stimulating gland normalization and promoting the thermal effects of OPE technology action,
- LIGHT MODULATION:
  - By emitting high-powered LED light, endogenous heat from both eyelids is triggered, helping to stabilize and increase the lipid tear film,
- SYNERGY:
  - Synergy between two medical technologies optimizes the treatment of dry eye-related pathologies,
- EYE-LIGHT PROTOCOL:
  - Treatment alleviates the symptoms of pathology in just a few hours,
  - Thanks to the synergy of the two technologies, Meibom's glands resume the production of lipids necessary for the proper functioning of the eye,
  - It is not painful, takes a few minutes and allows for an immediate return to normal activity,
  - You will notice an improvement in the function of the Meibom's glands just a few hours after the procedure,
- Modular screening screen:
  - Imaging technologies: 3D imaging, infrared imaging,
  - Screening time: usually completed within 5 minutes,
  - Integration: integration of Meiboscale and OSDI-6 tests for comprehensive evaluation,
  - Software: built using complex mathematical algorithms for accurate screening analysis,
  - User Interface: Intuitive and operator -friendly interface for ease of use,
- Lenovo Operator Laptop:
  - Laptop and 10.1-inch tablet in one, with detachable keyboard for enhanced capabilities,
  - Intel Atom x5-Z8350,
  - 64GB eMMC memory,
  - System: Windows 10,
- Camera placed on a mobile cart with 4 wheels including 2 with brake,
- Built-in wiring basket,
- Includes:
  - Power cord,
  - Operator glasses,
  - Patient glasses,
  - Eye lamp light module + head,
  - Treatment patient mask with light module,
  - Post-treatment patient mask with light module (also for dermatology services),
  - Operator's laptop/tablet,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

