



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**KRIOPOL R30 BRYZA II Cryotherapy unit (Ex-Display)**

- New (ex-display),
- Real photos of the product,
- Made in Poland, (2023),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- KRIOPOL R Bryza II device is designed for local cooling of tissue. The working medium is liquid nitrogen, thanks to which we obtain temperatures of -160° C at the nozzle outlet, which ensures effective and rapid stimulus stimulation,
- It is placed on a mobile frame, which allows it to be easily moved,
- It is equipped with a 30 dm3 liquid nitrogen tank,
- It has a flexible heat-insulated liquid nitrogen vapor transfer line, which makes it convenient and safe to apply coolness to any part of the patient's body,
- It has a segmented monochrome LED display panel,
- It has a four-stage regulation of blowing intensity, thanks to which we apply the right dose of coolness depending on the mass of the cooled area,
- It has additional blowing functions:
  - Release trigger points,
  - Two types of pulsed application,
  - Programmable treatment duration,
- The device is equipped with a constant measurement of the amount of nitrogen in the tank,
- The analgesic effect of low temperature allows for full kinesitherapy of the affected joints. The reduction of pain allows patients to perform exercises that they would not be able to perform in a painful state. The cryostimulation procedure is perfectly tolerated by patients, and the tissue cooling causes strong endogenous overheating of the joints until the baseline temperature is exceeded by about 3-4 °C. This effect lasts for about 3-4 hours,
- Treatment with low temperatures is followed by:
  - pain relief,
  - weakening of the activity of the inflammatory process,
  - reduction of muscle tension,
  - reduction of swelling,
  - clinical and functional improvement, which is an increase in the range of joint mobility,
  - increase in muscle strength,
  - alleviation of post-burn conditions,
  - reduction of injury healing time,
- Application:
  - Rehabilitation,
  - Sports medicine,
  - Wellness,
- Indications:
  - rheumatoid arthritis (RA),
  - ankylosing spondylitis(AS),
  - degenerative joint conditions,
  - cervical spine syndrome,
  - periarthritis,
  - painful shoulder syndrome,
  - discopathies,
  - Reiter's syndrome,
  - lichen erythematosus,
  - gunshot,
  - tendopathies,
  - myalgias,
  - tendon ruptures, peroneal tendonitis, tendonitis,
  - dislocations, sprains, fractures,
  - postoperative swelling,
  - burns,
  - long non-healing wounds,
- Included:
  - Power cord,
- Dimensions: 47 x 54 x 100 cm,
- Weight: 34.5 kg (without nitrogen),
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

