

PRAISTON Sp. z o.o. ul. Górowska 32 64-100 Leszno Polska Biuro: +48 65 527 01 67 Serwis:+48 782 844 000 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży Krzysztof Wybieralski +48 882 762 006

kw@praiston.pl



OLYMPUS WM-NP1 endoscope cart (Reconditioned)

Reconditioned (used), Technical condition: very good, Visual condition: very good, Real photos of the product, U.K. production (2013), Power supply: 230 V, Frequency: 50/60 Hz, Power: 1800 VA, Olympus WM-NP1 is a stable and sturdy cart suitable for carrying endoscopes and delicate medical equipment, It has: Upper tabletop with dimensions: 47 x 38 cm, Lower top with dimensions: 47.5 x 56 cm, 2 shelves with dimensions: 40 x 54 cm with manual height adjustment (every 1.5 - 2 cm) with Allen screws, Handles for easy movement of the cart, Holder for endoscopes with manual height adjustment, Monitor holder with height adjustment within 20 cm (6 points every 3.5 cm), Electrical installation, 12 power sockets for the medical equipment carried, Product features: Power switch on the front of the cart, Detachable cart back panel for easy wire management, Ability to hide wires in panels on the sides of the cart. Corners of shelves protected with protective pads and edges over wheels with rubber, Mobile on four wheels with Ø 12 cm, two of which have brakes, Shelf with Velcro fasteners designed for 2 bottles, Selectable menu in language:, Keypad on the control panel in the language:, Included: Power cord built-in transformer, MH-155 white balance cap, Overall dimensions: 67 x 66 166/200 cm, Weight: 86.5 KG, Cat 08. Has a current inspection and is ready for operation, Issued Technical Passport (Service Report) valid for a period of 12 months, Warranty: 6-months for domestic market (Poland), 3-months for the international market, Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge, Financing options (Poland only): Installments, Leasing, Loan,



AISTON









www.praiston.pl

2022