



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Magnetotherapy device ASTAR MAGNERIS (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production, (2013),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 250 VA,
ASTAR MAGNERIS

is a device designed for low-frequency magnetic field therapy. The device works with a modern two-part flat applicator. The magnetic field generated during therapy is concentrated between the parts of the applicator, which makes it possible to reduce its impact on the people operating the device, while increasing its use in the affected area,
The two-part flat applicator, which is included in the equipment, allows optimum use of the generated magnetic field by focusing it mainly on the lesion area. The advantage of such a solution is that the impact of the field on those operating the device and those in the vicinity is significantly reduced,
The controller and applicator are attached to a tripod equipped with casters, and the tripod is additionally equipped with a pneumatic actuator, balancing the weight of the application part, thus improving working comfort,

Application:
Delayed bone fusion,
Pseudarthritic joint,
Osteoporosis,
Joint degeneration,
Rheumatoid arthritis,
Ulcers and trophic changes of the lower extremities,
Bacterial infections of the skin and soft tissues,
Scarring,
Post-stroke conditions,
Migraine, vasomotor headache,
Cranial and peripheral nerve dysfunction,
Multiple sclerosis,
Infectious disease of the cornea,
Disappearance of the optic nerve,
Hypertension,
Coronary artery disease,
Cardiac arrhythmia,
Irritable colon,
Long-term pancreatitis,
Technical specifications:
Maximum induction: 0-2.5mT/0-25Gs,
Modulated mode parameters: pulse 1s/ interval 0.5-8s,
Frequency range: 2 - 60Hz,
Treatment timer: 1-30 min

Magnetic field:
Field shape: sine, triangle, rectangle, half-sine, half-triangle, half-rectangle,
Continuous and modulated emission,
Wide frequency range

Ergonomics:
Operating modes: program/manual,
Large, easy-to-read graphic display - clear operation in graphic mode,
Selection of disease units by disease name,
Treatment timer,
Two dedicated applicator slots,
Automatic recognition of the type of connected applicator,
Included:
Power cord,
CP-ID1 | CP-ID2 flat applicator,
Dimensions: 125 x 70 x 110 cm,
Weight: 25 KG,
Has up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

