



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Plasma Fission (New)

New,

A device to perform treatments using a plasma beam - it has many features that distinguish it from other devices. The treatment with this device is based on selective elimination of fat cells by means of focused ultrasound waves without disturbing surrounding tissues, blood vessels and nerves,

Indications:

Eye lift (non-invasive bluffing),

Lifting the lower part of the face,

Reduction of wrinkles in the eye and lip area (so-called smoker wrinkles),

Removal of skin lesions,

Main applications:

Removal of drooping eyelids,

Elimination of wrinkles over the upper and lower lips - the so-called "smoker's wrinkles",

Removal of micro-wrinkles in the area of eyes and ears (crow's feet),

Elimination of flaccid skin in the area of elbows and knees,

Removing all kinds of skin lesions:

Seborrheic warts,

hemangiomas,

fibroma,

...yellows,

the warts,

Correction of small, protruding scars,

Reducing skin discoloration (ostudes),

Removal of hyperkeratosis symptoms, i.e. excessive keratosis of the epidermis,

Used in acne treatment support therapy,

Included in the set:

MicroPlasma device,

Quick Charge USB 2.0 charger,

Aluminium case,

Treatment electrodes,

A tripod,

10 disposable electrodes,

Characteristics:

Easy to read battery level indicator,

Painless and safe for the patient,

The use of plasma does not cause inflammation, bleeding and necrosis of surrounding tissues, because its action is limited to the border of the basal layer of the dermis,

The plasma generated by the contact between the device's electrode and the patient's skin produces a thermal effect, which results in protein coagulation and tissue evaporation,

A great alternative to the scalpel,

MicroPlasma is fully wireless:

No passive electrode,

battery power supply,

Charging with USB charger in Quick Charge 2.0 standard,

The magnetic charging interface in the device protects the device against breakage and makes it easy to connect the charger quickly,

Energy saving - when not in use, the device will sleep and signal reduced power consumption, and after a certain period of time it will switch off,

The MicroPlasma tripod allows you to safely place your MicroPlasma between treatments and during charging,

The aluminium case will protect your MicroPlasma during transport,

Intuitive operation,

Colourful OLED display,

Ergonomic shape,

9-step power adjustment,

6-level security system,

Technical data:

Ergonomically designed to hold the device securely - the weight is only 135 grams,

Battery power supply,

Intuitive operation - clear, colourful OLED display

Nine-stage power control,

A plasma generator with automatic voltage regulation ensures constant beam power and safety during the treatment,

Extensive software and hardware security system,

The MagClick™ system for mounting the electrode in the head allows for quick replacement of the electrode,

Disposable sterile treatment electrodes with a diameter of 0.15 mm enable precise targeting of the plasma beam,

Built-in 8Wh battery allows you to perform treatments all day without the need for recharging,

He has issued a valid Technical Passport,

Warranty 24 months

