

PRAISTON Sp. z o.o. ul. Górowska 32 64-100 Leszno Polska

Biuro: +48 65 527 01 67 Serwis:+48 782 844 000 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży Krzysztof Wybieralski +48 882 762 006

kw@praiston.pl



## SORISA IMPRESS LITE MED lymph drainage device (Ex-Display)

New, exhibition, Spanish production, Power supply: 230 V, Frequency: 50/60 Hz,

Power consumption: 70 - 80 VA,

Presotherapy is a method widely used in the treatment of lymphatic and circulatory systems, in cosmetic medicine, physiotherapy and cosmetics,

Maximum safety thanks to the unique blood pressure monitor with APS (Adapted Pressure System), which automatically personalizes the parameters of each treatment,

Lower and upper limb mode with 10 filling modes and 5 adjustable pressure levels,

18 programmed upper and lower limb treatment programs,

Full selective swelling control (from one limb),

Color touch screen,

Clear and easy to use interface, Outlet pressure characteristics:

Without measuring blood pressure: 140 mmHg,

Performing a blood pressure measurement: 90% systolic pressure,

Maximum pressure regulation tolerance: ± 5 mmHg,

Blood pressure measurement:

Method: oscillometric Pulsation: 40 - 120,

Maximum systolic pressure: 180 mmHa. Minimum diastolic pressure: 40 mmHg,

treatments: Cellulite,

Tired legs and drainage, Circulation stimulation.

Lymphoedema and venous edema. Postoperative and traumatic edema, Cellulite and obesity prevention,

Postnatal care

General detoxification of the body,

The kit includes:

SORISA IMPRESS LITE MED lymph drainage device,

Jumpsuit for legs and feet with a cuff for measuring pressure,

Two upper limb cuffs,

Hip cuff, Power cord,

Dimensions: 400 x 350 x 230 mm,

Weight: 9.4 kg,

Has a valid Technical Passport issued,

Manufacturer's warranty unit:

6 months for the domestic market - Poland,


1 month on the international market,

Financing option (Poland only): Installments, Leasing, Loans,






















www.praiston.pl