



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Cassette for sterilization of small instruments PRESTIGE MEDICAL (Reconditioned)

Used, very good condition,
English production,
The cartridge is used to sterilize small instruments in selected PRESTIGE MEDICAL autoclaves,
Pull-out lid,
Holes placed on the flap and on the bottom of the case,
Internal dimensions: 5.5 x 17.5 x 3 cm,
External dimensions: 7.4 x 19.3 x 4.3 cm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan.

