



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



K-1 stethoscope (Ex-Display)

A new one, post-warrior,
P.R.C. production,
Flat, metal, unilateral stethoscope,
Most commonly used for pressure measurements using mechanical pressure gauges,
It has one 54 cm long drain,
Dimensions:
Total length: 78 cm,
Diaphragm diameter: 4.5 cm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan.

