



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CF-Q180AL video colonoscope (Reconditioned)

Used, very good condition,
Japanese production,
Direction of vision: Forward view,
Outer diameter of the tip: 12.8 mm,
External diameter: 13.2 mm,
Working channel diameter: 3.7 mm,
Working length: 1680 mm,
Deep down: 3 - 100 mm,
Maximum bending angle of the endoscope tip up and down: 180°/180°,
Maximum bending angle of the tip right and left: 160°/160°,
Maximum angle of field of vision: 170°,
Fully processor-compatible: CV-180/160/140,
It's complete:
Original suitcase,
Lid,
Valves,
User manual in English (PDF),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan,
OVERVIEW PICTURES OF THE DEVICE

