



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Orthopaedic flange NECLOC NL-200E (Reconditioned)

Used, good technical condition,

U.S.A. production,

Size: Pediatric (red),

NecLoc is the world's most widely used two-piece orthopaedic collar, used in the extraction of accident victims,

It is a fast and safe system to immobilize cervical vertebrae,

The basic structure is made of low-density polyethylene, whose flexibility makes it easy to use and rigidity optimises the ability to immobilise,

The patient's comfort is ensured by the use of an internal padding made of expandable, closed-cell, high-density plastisol,

The anterior hole allows the tracheotomy to be performed,

Transparent for X-rays,

Easy to keep clean,

Each size is marked with a different color to allow quick identification of the correct flange,

Cat 02,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan.

