



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Pinhole cladding (blue) (Ex-Display)

A new, post-exposure one,
P.R.C. production,

Pinhole cladding is used in ophthalmological surgeries to examine the "pinhole" acuity,

The pinhole test is a very quick way to get more information about a person's eye condition,

The cladding is made of plastic and consists of a handle and a plate covering one eye. The other eye is covered by a tiltable flap with tiny holes. There is a groove on the flap for the patient's nose,

Dimensions:

Handle length: 13 cm,

Total length: 20.3 cm,

A plate that covers the eyes: 15,2 x 6,6 cm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

