



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Videolaryngoscope Pocket for children (PVL-I-01) (New)

New,

PRC production,

The Pocket Videolaryngoscope is designed for everyday use in conventional and complex intubation procedures,

The pocket Videolaryngoscope is designed to provide physicians with a consistent, clear picture of the patient's airway in real time and to help them achieve faster and safer endotracheal intubation,

The model (PVL-I-01) is designed to work with children and infants,

Technical data:

Monitor type: Bi-directional rotating / fully assembled,

Size / ratio of the monitor: 3"/4:3,

Viewing angle: 85°,

Rotate the monitor: Two-way left/right rotation $\geq 300^\circ$; front/rear $\geq 180^\circ$,

Camera pixels: 4500000P CMOS,

Frame rate: 30 fps,

Field of vision (FOV): $65 \pm 10^\circ$,

Depth of field (DOF): 10-80mm,

Internal power supply / battery: 18500 / 3.7 V DC / 2040 mAh replaceable / battery,

Battery life of a fully charged battery: ≥ 200 min,

Water resistance: IPX7 for VL mount; IPX4 for monitor,

Internal memory: 16G,

Operating temperature: 0°C-50°C,

Moisture: 30%-80%,

Storage temperature: -20°C-50°C,

Storage humidity: 93%,

Emergency power supply: NO,

Data recovery: Waterproof USB port,

Application: For paediatric patients,

Specialty: Pocket Portable / patented spoons,

Type of blades: Disposable: MD-0, MD-1, MD-1S,

Product color: Black,

He has a current Technical Passport issued,

Warranty 12 months,

Possibility of financing (Poland only): Installments, Leasing, Loan.

