



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Dental mirror HAHNENKRATT SEplus size 8 (Ex-Display)

New, post-exposure,
Germany production,

SEplus dental mirrors are made of surgical steel (stainless) and high quality laser cut glass, which ensures comfort and professionalism,
They guarantee perfect reflection, do not distort the image,

The mirror is suitable for sterilization in an autoclave or ultrasonic cleaner,

After repeated sterilization, the mirror is still usable and does not lose its initial parameters,

Mirror size: 4,

Dimensions:

Length of the working tip: 28 mm,

Mirror diameter: 29 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

