



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Dental mirror PDI size 4 (Ex-Display)**

New, post-exposure,  
U.S.A. production,

PDI dental mirrors are made of surgical steel (stainless) and high quality laser cut glass, which ensures comfort and professionalism,  
They guarantee perfect reflection, do not distort the image,

The mirror is suitable for autoclave sterilization at 134°C and 18 minutes sterilization time,

The mirror can also be sterilized with dry hot air at 160°C and a sterilization time of 2 hours,

After repeated sterilization, the mirror is still usable and does not lose its initial parameters,

Mirror size: 4,

Dimensions:

Length of the working tip: 30 mm,

Mirror diameter: 22 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

