



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RTG REGO cassette 20 x 40 cm (Reconditioned)

Used, very good technical condition,

Germany production,

REGO cassette is a plastic-metal cassette with a magnetic system allowing to obtain a perfectly sharp image. Magnetic system allows to obtain uniform and accurate film pressure to the screens on the whole surface. This makes the REGO cassette an ideal solution, e.g. for exposing the whole spine or lower limbs,

Size of membrane: 20 x 40 cm,

General dimensions: 22,8 x 42,8 x 1,5 cm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

