



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Hospital laundry cart (folding) (Reconditioned)**

Used, very good technical condition,  
Polish production,

Laundry trucks are designed to transport dirty laundry, used in hospital wards, sanatoriums or other medical facilities,

The cart is made of acid resistant stainless steel tubes,

It has a folding frame equipped with four casters with rubber bumpers,

Special hooks prevent the bag from slipping off the rim,

No bag included,

Dimensions of the unfolded cart: 60 x 61 x 100 cm,

Dimensions after folding: 20 x 60 x 115 cm,

Weight: 6.8 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

