



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LUXO T88/87 diagnostic lamp (Reconditioned)

Used, very good technical condition,
Italy production,
Power supply: 230 V,
Frequency: 50/60 Hz,
G23 bulb with 11 W power,
Single focal,
Fixed to the tabletop,
Bowl adjustable in many planes,
Adjustable articulated arm and rotated on the handle,
Lamp switch placed on the canopy,
The lifetime of the bulb: 10000 hours,
Color rendering: 80-89 Ra,
Luminous flux: 900 lumens,
Color temperature: 3000 K,
Includes mounting bracket and power supply unit,
Canopy dimensions: 9 x 33 x 5 cm,
Length of arms after straightening out: 80 cm,
Weight of the set: 1.6 kg,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan.

