



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Bone screws DANEK (Ex-Display)

New, post-exposure,

U.S.A production,

Two bone screws included:

Screw DANEK 823-110:

Thread diameter: 6.5 mm,

Diameter of the head: 8 mm,

Thread length: 25 mm,

Total length: 40 mm,

Screw DANEK 823-112:

Thread diameter: 6.5 mm,

Diameter of the head: 8 mm,

Thread length: 38 mm,

Total length: 50 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

