



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical screw kit (Reconditioned)

Used, very good technical condition,

Included:

Screw 808-945 41830:

Thread length: 45 mm,

Thread diameter: 7.5 mm,

Total length: 60 mm,

Screw 808-945 87118:

Thread length: 45 mm,

Thread diameter: 7.5 mm,

Total length: 60 mm,

Screw 808-945 82952:

Thread length: 45 mm,

Thread diameter: 7.5 mm,

Total length: 60 mm,

Screw 808-945 60492:

Thread length: 45 mm,

Thread diameter: 7.5 mm,

Total length: 60 mm,

Screw X0194-005 88152:

Thread length: 35 mm,

Thread diameter: 8.5 mm,

Total length: 50 mm,

Screw 5:

Thread length: 33 mm,

Thread diameter: 5 mm,

Total length: 49 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

