



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Set of long piston blades STRYKER (Reconditioned)

Used, very good technical condition ,
U.S.A production,
Model: 277-96-300,

The set includes four blades with dimensions:

Blade length: 77.5 mm,

Blade width: 6.3 mm,

Blade thickness: 0,64 mm,

Cutting thickness: 1,04 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

