



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ORION PLATE HOLDING PIN DRIVER Surgical Screwdriver (Reconditioned)

Used, very good technical condition ,
Model: 871-655,
Orthopaedic instrument made of stainless steel,
It has a wooden handle with a screw sleeve for surgical and veterinary use,
Dimensions:
Shank diameter: 5,8 mm,
Shank length: 15 cm,
Screw length: 10 mm,
Total length of the screwdriver: 26,5 cm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan.

