



### **NESLAB RTE-101 cooling water bath (Reconditioned)**

Used, very good technical condition,

U.S.A production,

Power supply: 230 V,

Frequency: 50/60 Hz,

The water bath with cooling offers efficient cooling and heating over a wide temperature range and very good water circulation,

Temperature range: -12° to +130°C,

Pump capacity: max. 11.5 l/min,

Tank capacity: 4.5 liters,

Tank dimensions: 12,7 x 12,7 x 12,7 cm,

General dimensions: 57,8 x 25,7 x 36,8 cm,

Weight: 29.5 kg,

User manual in English (PDF),

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

