



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER M-1 ambulance stretcher (Reconditioned)

Used, very good technical condition,
U.S.A. production (2014r.),
Maximum load: 227 kg,
It has a lightweight, durable construction made of light alloys,
Conveyor frame made of rectangular profiles (increased resistance to extreme overloads),
Adapted for resuscitation hard ground,
Multi-level stretcher with independent front and rear adjustment (6 levels),
Possibility to set the conveyor in a shock-absorbing position,
Trendelenburg and Fowler position,
A set of safety belts - harness, cross and paediatric belts,
Retractable side rails,
Pull-out handles for carrying,
A system of fast and secure connection with stretchers, allowing the stretcher to be fastened forwards and backwards,
Possibility of leading the stretcher sideways by one person from any side of the stretcher,
Automatic chassis folding system for easy loading into the ambulance,
4 castors - all 360° swivel castors with the possibility of locking 2 castors into fixed position,
Brakes on 2 wheels,
Double protection against accidental folding of the conveyor chassis,
Dimensions: 55 x 198 x 110 cm,
Weight: 66,5 kg,
The user manual in Polish, English and German (PDF),
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

