



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Dental ball pad FM size 2 (Reconditioned)

Used, very good technical condition ,
Polish production,
Made of medical stainless steel,
Classic double-sided, double-ended ball pad,
The tool has an anti-slip handle, which guarantees good tool guidance and good grip,
It is used to develop medium size fillings,
Tool length: 170 mm,
Tips with diameter:
1,57 mm,
1,96 mm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

