



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TESON 5.0/7.5 MHz endovaginal head (Reconditioned)

Used, very good technical condition,

Polish production,

Type: endovaginal,

Frequency: 5.0/7.5 MHz,

Dimensions:

Head (handle): 3,7 x 4,5 x 27,5 cm,

Cable length: 220 cm,

Included in the set of suitcase,

Cat 02,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

