



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Linear head ECHOSON LT-75 (Reconditioned)**

Used, very good technical condition,

Polish production,

The name of the probe: LT-75,

Type: linear,

Frequency: 7.5 MHz,

Probe compatible with ECHOSON cameras (Albit, Spinel, Spinel II),

It provides very good imaging to a depth of about 10 cm,

Ideal for veterinary, endocrinology purposes, it allows for differential diagnosis of very small morphological (focal) lesions, as well as a very detailed diagnosis of subcutaneous changes,

Dimensions:

Head (handle): 7.9 x 3 x 12.3 cm,

Cable length: 170 cm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

