



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Linear head ECHOSON LT-75 (Reconditioned)

Used, very good technical condition,

Polish production,

The name of the probe: LT-75,

Type: linear,

Frequency: 7.5 MHz,

Probe compatible with ECHOSON cameras (Albit, Spinel, Spinel II),

It provides very good imaging to a depth of about 10 cm,

Ideal for veterinary, endocrinology purposes, it allows for differential diagnosis of very small morphological (focal) lesions, as well as a very detailed diagnosis of subcutaneous changes,

Dimensions:

Head (handle): 7.9 x 3 x 12.3 cm,

Cable length: 170 cm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

