



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Dental canal wound hook (Reconditioned)**

Used, very good technical condition,  
Germany production,  
High quality surgical hook made of German stainless steel,  
The high quality of workmanship and the materials used guarantee functionality, durability and high resistance to corrosion under conditions of use, disinfection and sterilization,  
The hook can be used in general and dental surgery procedures,  
The product is designed for repeated use after previous washing, disinfection and sterilization,  
Total length: 270 mm,  
Handle width: 14 mm,  
The width of the hook tip: 8 mm,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Installments, Leasing, Loan

