



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video - colonoscope OLYMPUS Evis Exera CF-Q145L (Reconditioned)

Used, very good technical condition,
Japanese production,
Direction of vision: Forward view,
Maximum field of view angle: 140o,
External diameter: 12.8 mm,
Working channel diameter: 3.7 mm,
Working length: 1680 mm,
Total length: 2000 mm,
Depth of focus: 3-100 mm,
Max. The bending angle of the endoscope tip up and down: 180o/180o
Max. Bending angle of the tip right and left: 160o/160o,
Fully compatible with CV-180/160/140 processors,
Complete:
Original suitcase with dimensions: 60 x 43 x 13,5 cm,
Valves,
Lid,
It has a current overview, and a Technical Passport,
6-month warranty,
OVERVIEW PICTURES OF THE DEVICE

