



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Round tip tweezers BD 765 (Reconditioned)

Used, very good technical condition,
Germany production,

High quality precision microsurgical tweezers made of German non-oxidizing stainless steel (INOX),
The high quality of workmanship and the materials used guarantee functionality, durability and high resistance to corrosion under conditions of use, disinfection and sterilization,

Tweezers are mainly used for grabbing and taking small pieces of tissue,

The anatomical notches on the round beak enable sensitive tissues to be gripped without risk of damage,

The product is designed to be reusable after previous washing, disinfection and sterilization,

Dimensions:

Total length: 200 mm,

Bow diameter: 4,5 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

