



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ESCHENBACH MAX TV magnifying glasses (Reconditioned)

- Used, very good technical condition,
- Germany production,
- MaxTV magnifying glasses are designed mainly for visually impaired people with visual acuity of 0.2 - 0.4,
- Magnification: 2.1x,
- Range of pupil spacing: 60 - 68 mm,
- The optics can be individually adjusted in focus: ± 3 dioptries,
- The earmuffs can be adjusted without warming up,
- Original case included,
- Dimensions:
 - Length: 18.5 cm,
 - Width: 14.5 cm,
 - Height: 3 cm,
- Weight: 49 g,
- Version A,
- He has a current Technical Passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 1 month for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan

