



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



ESCHENBACH MAX TV magnifying glasses (Reconditioned)

Used, very good technical condition,
Germany production,
MaxTV magnifying glasses are designed mainly for visually impaired people with visual acuity of 0.2 - 0.4,
Magnification: 2.1x,
Range of pupil spacing: 60 - 68 mm,
The optics can be individually adjusted in focus: ± 3 dioptries,
The earmuffs can be adjusted without warming up,
Original case included,
Dimensions:
Length: 18,5 cm,
Width: 14,5 cm,
Height: 3 cm,
Weight: 49 g,
Version B,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

