



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Sectoral cardiology head HP 21242A (Reconditioned)

Used, very good technical condition,

U.S.A. production,

The name of the probe: 21242A,

Type: sectorial (for cardiac applications in adults),

Frequency: 3.5 / 2.7 MHz,

The probe is compatible with HP Sonos 1000, 1500 and 2000 series ultrasonic systems,

Dimensions:

Head (handle): 4,3 x 3 x 8,7 cm,

Cable length: 180 cm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

