



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Single-sided straight alveolar spoon (medium) (Reconditioned)**

Used, very good technical condition,

German production,

Single-sided straight alveolar spoon made of top quality stainless steel. It is used during surgical procedures in which there is a need to remove tooth splinters, ventricular septum, granules or granules,

The surgical teaspoon is finished on one side, unlike a bent tip teaspoon, it is designed for the jaw and not for the mandible,

The tool is profiled and lies perfectly in the hand, thus increasing safety during treatments. Thanks to the grooved surface it lies securely in the hand and does not slip out of it,

Dimensions:

Overall length: approx. 150 mm,

Spoon: 3,9 x 7 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

