



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Laboratory refrigerator LABCOLD RLDF0510A (Reconditioned)**

Used, very good technical condition,

English production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 95 W,

Professional laboratory fridge with a capacity of 150 litres powered by a unique IntelliCold controller, designed from scratch for safe storage of medicines in the cold chain,

Equipped with IntelliCold controller,

Automatic temperature recording 24/7,

Power failure alarm,

Too high a temperature alarm,

Too low a temperature alarm,

Open door alarm,

Automatic defrosting with switch-off possibility,

Fan assisted cooling,

Environmentally friendly refrigerant,

Internal LED light,

Capacity: 150 litres,

Temperature range: +2 to +8°C,

Operating temperature: +5°C,

Three shelves with height adjustment,

Locks the door,

External dimensions: 835 x 600 x 600 mm,

Internal dimensions: 690 x 480 x 410 mm,

Weight: 51 kg,

The user manual in English and German (PDF),

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

