



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Dental scaler (sickle) folded (small) (Reconditioned)

Used, very good technical condition,
German production,
Designed for manual removal of deposits, plaque and tartar,
The tool is profiled and lies perfectly in the hand, thus increasing safety during treatments. Thanks to the grooved surface it lies securely in the hand and does not slip out of it,
Made entirely of surgical steel,
Possibility of sterilization in an autoclave,
Dimensions:
Overall length: 153 mm,
The length of the tip: 9 mm,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

