



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Hemingway two-sided (small) Hemingway spoon Version B (Reconditioned)

Used, very good technical condition,

Polish production,

Made of medical stainless steel,

Hemingway spoon is used to remove the grain from the alveolus. It is also used for exfoliating inflammatory lesions and removing broken parts of teeth and the interroot septum,

The bone cot consists of an anti-slip handle and two working parts. The right and left working part guarantees easy handling of the tool,

Tool length: 170 mm,

Spoon width: 2 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

