



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



One-sided FM dental curette (Reconditioned)

Used, very good technical condition,

Polish production,

Universal dental curette for preventive treatments in dentistry. The instrument is used to remove subgingival dental deposits (tartar removal).

It is also used in the final phase of removing supragingival deposits and smoothing the root surface using a mechanical method. This type of dental curettes complement the work of a dental scaler,

Made of stainless steel,

It is suitable for autoclave sterilisation,

Dimensions:

Overall length: approx. 140 mm,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

