



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Smear test / tests Panbio COVID-19 Ag ABBOTT 25 pcs / pack (New)**

New,  
German production,  
The Panbio COVID-19 rapid antigen test is a spot test to help detect an acute SARS-CoV-2 infection,  
The test is carried out on material taken from a nasopharyngeal swab,  
The result of the examination is available after only 15 to 20 minutes, so that the examination can be carried out directly on the spot in the doctor's office, hospital or clinic,  
The Panbio COVID-19 rapid antigen test is on the BfArM list and is therefore reimbursed in accordance with the "Regulation on the right to be tested in connection with the direct detection of the SARS-CoV-2 coronavirus pathogen",  
The packaging is included:  
Test cassettes - 25 pieces,  
Buffer solution (9 ml) - 1 piece,  
Tubes and lids for extraction - 25 pieces each,  
Smear with positive control - 1 piece,  
Smear with negative control - 1 piece,  
Sterile swabs to take samples from nasopharynx - 25 pieces,  
Extraction tube rack - 1 piece,  
Shortened guide - 1 piece,  
Multilingual instruction manual - 1 piece,  
Storage temperature: 2° - 30° C,  
Sensitivity: 93.3%,  
Specificity: 99.4%,

