



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Cryosurgery apparatus - portable Freezpen CL-FP-35 (New)**

New,  
Belgian production,  
Freezpen is manufactured by Clinic6,  
It is an innovative device increasing the overall improvement of treatment quality through direct application of Nitric Oxide with a temperature of  $-89^{\circ}\text{C}$ , penetrates in 1 mm for 5 seconds to a maximum depth of 5 mm,  
Ease of use is an incentive for family doctors, aesthetic medicine and paediatricians,  
Freezpen offers stress-free treatment with quickly visible improvement effects,  
The high precision of Freezpen reduces damage to healthy skin, minimising discomfort for patients,  
Effective result obtained thanks to the first treatments, thanks to the high pressure of the Freezpen preparation reaching deeper penetration,  
The camera has a calming pen shape, the advantage is optimal safety with minimal discomfort and no side effects,  
Very economical with cheap cryogenic capsules that can treat up to 12 patients,  
The camera uses cartridges weighing 35 grams,  
The device is an optimal solution for treatment:  
papillae,  
stabbing,  
hemangiomas,  
skin markers,  
corneas and all other mild skin disorders.  
The apparatus is included:  
3 applicators:  
1 applicator with green coding (for treatment of lesions with diameter  $\varnothing$  4-8 mm)  
1 applicator with orange - 60 mm (for treatment of lesions with diameter  $\varnothing$  1-4 mm),  
1 applicator with white coding -60 mm (for treatment of lesions with diameter  $\varnothing$  2-6 mm),  
Instructions for use,  
Tweezers for filter insertion.  
It has a current Technical Passport issued,  
Warranty 12 months

