



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Cable PIGTAIL OLYMPUS MAJ-1430 (Ex-Display)

The new post-exposure, perfect technical condition,
Japanese production,

It provides a connection between the video processor and the endoscope,

MAJ-1430 is designed for use with OLYMPUS CV-180 and CV-190 video processor and endoscopes series 180 and 190,
Spiral, flexible cable,

Weight: 0.8 g,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

