



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Pulse oximeter per finger Viatom OXYSMART PC-60F (New)**

New,  
P.R.C. production,  
Power supply: 3.0 VDC (2 x AAA),  
Menu in English,  
Pulse oximeter PC-60F per finger, is designed to measure heart rate and functional oxygen saturation by the patient's finger. It is used for random control of SpO2 and heart rate of adults and children in homes and clinics,  
Lightweight, compact and elegant design provides a stylish design,  
Quick and easy to use, one button operation,  
Finger-operated,  
Resistant to splashes and falls,  
Data, including SpO2 and heart rate, can be monitored in 5 seconds,  
OLED colour display with displayed functions:  
Oxygen saturation,  
Pulse,  
A pulse wave,  
Pulse bar graph,  
Battery indicator,  
Functions:  
Automatic on/off,  
4-way display,  
Sound and visual alarms,  
Alarm for incorrect readings,  
Low charge indicator baterii  
Point control and continuous measurement mode,  
Heart rate analysis,  
Memory for up to 12 SpO2 data storage groups,  
Specifications:  
Pulse frequency range: 30 to 240 beats per 1/8,  
Oxygen level range: 70%-100%,  
Perfusion (pi) tange: 0%~20%,  
Oxygen transmitter: double wavelength led sensor,  
Measurement mode: 30-second spot check and continuous mode,  
Default alarm limit Spo2 low limit: 90%,  
Default limit for heart rate alarm:  
High limit: 120 bpm,  
Low limit: 50 bpm,  
Included:  
Pulse oximeter,  
Cover,  
Lanyard to hang,  
The user manual in English,  
1.5 V (AAA) batteries - 2 pcs,  
Dimensions: 34 x 30 x 56 mm,  
Weight: 52 g,  
He has a current Technical Passport issued,  
Warranty 12 months

