



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Surgical ophthalmic magnifier PENTAX PRISMVUE NF2 T Series S (3.0x) (New)

New,

Japanese production,

The S-series magnifier is a high-end Bright model,

PRISMVUE magnifier - is a prismatic magnifier equipped with unique PENTAX deflection prisms which provide the surgeon with a bright, enlarged view of the surgical field, while maintaining a straight posture, preventing the need to tilt the neck and back forward

The vertical forward-facing posture during operation prevents the operating lights from interfering with the user's head and always provides a bright field of vision.

Advantages of the magnifying glass:

Unique optical design,

A clear and broad field of vision,

The field of vision is 20-30% larger than that achieved with conventional products,

Reduces eye fatigue and neck tension,

Easy switching between enlarged and real views,

Replaceable noses,

Technical data:

Magnification: 3.0 x (105g),

Lens configuration:

3 Lens,

3 Spectacle lenses,

Total 6,

Working distance (WD): on request (35 cm-70 cm),

Frame: Nylon

Available frame colours:

Red,

Blue,

Black,

Silver,

He has a current Technical Passport issued,

24 month guarantee,

Possibility of financing (Poland only): Instalments, Leasing, Loan.

