



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Surgical ophthalmic magnifier PENTAX PRISMVUE NF2 L Series C (3.5x) (New)**

New,

Japanese production,

The C-series magnifier is a high-end ultra-light model,

PRISMVUE magnifier - is a prismatic magnifier equipped with unique PENTAX deflection prisms which provide the surgeon with a bright, enlarged view of the surgical field, while maintaining a straight posture, preventing the need to tilt the neck and back forward

The vertical forward-facing posture during operation prevents the operating lights from interfering with the user's head and always provides a bright field of vision.

Advantages of the magnifying glass:

Unique optical design,

A clear and broad field of vision,

The field of vision is 20-30% larger than that achieved with conventional products,

Reduces eye fatigue and neck tension,

Easy switching between enlarged and real views,

Safety glasses,

Technical data:

Magnification: 3.5 x (93g),

Lens configuration:

2 Lens,

3 Spectacle lenses,

Total 5,

Working distance (WD): on request (35 cm-70 cm),

Frame: Nylon

Available frame colours:

Red,

Blue,

Black,

Silver,

He has a current Technical Passport issued,

24 month guarantee,

Possibility of financing (Poland only): Instalments, Leasing, Loan.

