



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEPU MEDICAL SARS-CoV-2 Certified Antigen Test(s) (25 pcs/pack) (New)

New,
manufactured by P.R.C.,
It has an approved CE mark,
SARS-CoV-2 antigen rapid test kit (colloidal gold immunochromatography),
Product designed for qualitative detection of the presence of SARS-CoV-2 antigens in clinical samples (nasal swab),
Card-like format,
The test card is based on a specific antibody-antigen reaction. Contains colloidal gold labelled with SARS-CoV-2 N monoclonal antibody,
Simple to use,
Convenient, no equipment required,
Rapid, results in 15 minutes,
Stable, with high accuracy,
Storage temperature: 4-30 °C,
Specificity: 99.26 %,
Package includes:
Certificate,
Test cards - 25 pieces,
Swab sticks - 25 pieces,
Bottles with sample preparation solution (2.5 ml) - 2 pieces,
Instructions for use in Polish, English, Russian, Spanish and Portuguese

