



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybierański
+48 882 762 006

kw@praiston.pl



CLUNGENE Certified SARS-CoV-2 Antigen Test / Cassette Test (25 pcs / pack) (New)

New,
manufactured by P.R.C.,
It has an approved CE mark,
Lungene SARS - CoV-2 is a COVID-19 Antigen rapid test cassette for use by trained clinical laboratory personnel,
It is a lateral flow immunoassay for the qualitative detection of SARS-CoV-2 antigens in nasopharyngeal and oropharyngeal swabs,
Test result available after 15 minutes,
Storage temperature: 4-30 °C,
Sensitivity: 98.5 %,
Package includes:
Certificate,
Test cassettes - 25 pieces, (each cassette comes in a package with desiccant),
Extraction reagent ampoules (0.3ml) - 25 pieces,
Sterile disposable sticks - 25 pieces,
Test tubes - 25 pieces,
Test tube rack - 1 item,
Instructions for use in Polish, English and German

