



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Light source RICHARD WOLF RIWO LP 4008 (Reconditioned)

Used, very good technical condition,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Light source - two separate halogen bulbs: 24 V / 150 W,
Light bulb selection switch,
Two-step brightness control,
Built-in colour filters:
Red,
Orange
Yellow,
Equipped with a carrying handle,
Included:
Optical fibre R. Wolf 8061.35 length: 1,8m,
Light guide handle with holder for tube R. Wolf 8836.90,
Light guide tube R. Wolf 8836.08, length: 8 cm, diameter: 2.4 cm,
R. Wolf tube 8836.04 with a length of 11 cm and a diameter of 1.8 cm,
Tubus R. Wolf 8836.01 with a length of 13.5 cm and a diameter of 2.2 cm,
Power cable,
Dimensions: 24.5 x 16.5 x 10/13 cm,
Weight of set: 6.5 kg,
Cat 03,
Has a current Technical Passport issued,
Warranty:

6 months for the domestic market - Poland,
1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

If you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

