



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **ALARIS GATEWAY workstation (basic version) (Reconditioned)**

Used, very good technical condition,  
Swiss made, (2014),

Power supply: 230 V or via built-in NiMH 9.6 V/ 2.4 Ah battery,

Charging time:

When no pumps are connected: up to 8 hours to a charge level of 95%, depending on the battery consumption,

When pumps are connected: up to 16 hours to a charge level of 95%, depending on battery life and connection of 11 pumps ,

Operating time: 20 minutes,

Frequency: 50/60 Hz,

Power: 460 VA,

Alaris Gateway Workstations are designed as a modular system to provide communication between an Alaris infusion pump and any patient data management system (PDMS), patient monitoring system (PM), hospital information system (HIS) or clinical information system (CIS) requiring access to infusion data stored on the pump,

No pumps included - can be retrofitted with all ALARIS pumps from our range at additional cost,

Alaris stations come in three different configuration options, this one is the basic version with organisational functions and central AC power function. It does not support communication with the PDMS system. It can be used as an auxiliary workstation to increase the number of pumps needed at the bedside,

It has space for 7 ALARIS infusion (syringe) pumps and for 3 volumetric (volume) infusion pumps,

Eliminates bundles of power cables by using a single power cord,

Audible and visual alarms,

IPX1 - protection against water droplets falling vertically from above,

Connectors:

RS232,

Nurse call

Ethernet networks,

Auxiliary socket,

Included:

Two mounting brackets,

Instruction manual in Polish (PDF),

Power cable,

Dimensions: 49 x 125 x 22 cm,

Weight: 15 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan

